

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY TEMPLATE**

A. 510(k) Number:

k073071

B. Purpose for Submission:

New device

C. Measurand:

Genotype of Cytochrome P450 2C9 (CYP450 2C9) and Vitamin K epoxide reductase complex subunit 1 (VKORC1)

D. Type of Test:

Genotyping Realtime PCR

E. Applicant:

TrimGen Corporation

F. Proprietary and Established Names:

eQ-PCR™ LC Warfarin Genotyping kit

G. Regulatory Information:

Product Code	Classification	Regulation Section	Panel
Cytochrome P450 2C9 (CYP450 2C9) Drug Metabolizing Enzyme Genotyping System (ODW)	Class II	21 CFR 862.3360 Drug metabolizing enzyme genotyping system	Toxicology (91)
Product Code	Classification	Regulation Section	Panel
INSTRUMENTATION FOR CLINICAL MULTIPLEX TEST SYSTEMS (NSU)	Class II	21 CFR 862.2570 - INSTRUMENTATION FOR CLINICAL MULTIPLEX TEST SYSTEMS	CLINICAL CHEMISTRY
VITAMIN K EPOXIDE REDUCTASE COMPLEX SUBUNIT ONE (VKORC1) GENOTYPING SYSTEM (ODV)	Class II	21 CFR 864.7750 - PROTHROMBIN TIME TEST	HEMATOLOGY

H. Intended Use:

1. Intended use(s):

Refer to Indications for use below.

2. Indication(s) for use:

The eQ-PCR™ LC Warfarin Genotyping kit is an in vitro diagnostic test for the detection and genotyping of two single nucleotide polymorphisms (SNP) in the cytochrome P450 enzyme gene CYP2C9 known as CYP2C9*2 (C430T) and CYP2C9*3 (A1075C), and a SNP in the vitamin K epoxide reductase complex 1 gene VKORC1, known as VKORC1 (-1639G>A) obtained from EDTA – anticoagulated whole blood samples. The eQ-PCR LC Warfarin Genotyping kit is used as an aid in the identification of patients at risk for increased warfarin sensitivity. It is a qualitative assay for use in clinical laboratories upon prescription by the attending physician.

3. Special conditions for use statement(s):

For Prescription use only

4. Special instrument requirements:

For use with LightCycler® Real-Time PCR System instrument model 1.2.

I. Device Description:

eQ-PCRTM LC Warfarin Genotyping Kit reagents

The reagents in the eQ-PCR™ LC Warfarin Genotyping Kit are designed to be used specifically with the Roche LightCycler®. They are provided for 32 tests per kit and include:

- The **PCR Mix** is a pre-mixed reagent for PCR amplification including a DNA polymerase, PCR Buffer and dNTPs.
- The **LC1** is the Probe Mix 1 containing primer and probes designed to detect the 2C9*3 (A1075C) SNP in the CYP2C9 gene and the (-1639 G>A) SNP in the VKORC1 gene.
- The **LC2** is the Probe Mix 2 containing primers and probes designed to detect the 2C9*2 (C430T) SNP in the CYP2C9 gene.
- The **NF Water** is nuclease-free water.

J. Substantial Equivalence Information:

1. Predicate device name(s):

Verigene® Warfarin Metabolism Nucleic Acid Test

2. Predicate 510(k) number(s):

k070804

3. Comparison with predicate:

The following table summarizes the similarities and differences between the new and predicate devices.

Similarities:

Characteristics	Nanosphere Warfarin Assay	eQ-PCR™ LC Warfarin Genotyping Kit
	Predicate device (K070804)	Proposed device
Instrument	Verigene system with Verigene cartridge	Roche LightCycler with glass capillary tubes
Method to determine specificity	Hybridization	PCR Hybridization Melting Curve Analysis
Assay conditions	A closed system	A closed tube assay
Procedures	Reaction and SNP analysis require separate operations (transfer cartridge to reader)	PCR and SNP analysis in same tube (no transfer steps)

Differences:

Gene	CYP2C9 and VKORC1	CYP2C9 and VKORC1
Specimen type	Purified DNA from human whole blood samples	Purified DNA from human whole blood samples

K. Standard/Guidance Document Referenced (if applicable):

Drug Metabolizing Enzyme Genotyping System – Class II Special Controls Guidance Document

L. Test Principle:

The eQ-PCR™ LC Warfarin Genotyping Assay uses real time PCR technology to amplify and detect the target DNA. The extracted DNA sample is mixed with a PCR Mix, and an eQ-PCR™ specific Probe Mix reagent containing specific primers and fluorescent labeled probes for the CYP2C9 and/or VKORC1 gene polymorphisms. Amplification and detection are then performed in the Roche Diagnostics LightCycler® Real-Time PCR System instrument model 1.2 using conditions defined in the specific eQ-PCR™ LC Warfarin Genotyping Kit labeling. After the PCR reaction is completed, the LightCycler System automatically proceeds to the melting curve-based detection method. The signal detection is based on fluorescence resonance energy transfer (FRET) technology. Two probes are designed for each SNP

and labeled with different fluorophores. The anchor probe selectively hybridizes to the sequence flanking the SNP site, and the SNP detection probe recognizes the SNP and hybridizes to the sequence containing the SNP site. During FRET, the fluorophore on the anchor probe is excited by the light source of the LightCycler and the excitation energy is transferred to the fluorophore on the SNP detection probe. The emitted fluorescence is then measured at the respective wavelength.

The SNP detection is performed by melting curve analysis, a method that discriminates different allelic forms of DNA by melting temperature. The SNP detection probe is designed to have a unique melting temperature to the wild type and mutant SNP. With change in temperature, the probe dissociates from the target sequence at a specific melting temperature depending on the allelic composition of SNP. The SNP type is then detected based on the melting curve specific to each SNP.

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. Precision/Reproducibility:

To evaluate the within-laboratory, between operator, and between laboratory reproducibility of the eQ-PCR™ LC Warfarin Genotyping Kit on the Roche LightCycler®, three sites (two external and the sponsor's laboratory) and five operators participated. A panel of 31 samples was tested at each of the three sites on each of five days. At the two external sites, two operators per site performed the testing on alternate days. The panel consisted of whole blood (WB) samples collected in EDTA from 20 donors, and 9 genomic DNA samples (cctl below). These samples represented all the genotypes tested by this device. Each site used a different DNA extraction method for isolating DNA from whole blood samples. Each site then tested the resulting extracted nucleic acids using the eQ-PCR™ LC Warfarin Genotyping Kit on the Roche LightCycler®. Genotypes for each sample were determined by bi-directional sequencing and compared with the results from the device.

The composition of the reproducibility panel with the call rate for all three sites combined for each sample is described in Table below. There were 6 no calls, but no incorrect calls reported. The entire study produced a correct call rate of 98.7%.

Genotypic composition and percent call rate per sample:

Sample ID	Genotype			# samples tested	# genotype calls made	# correct calls	# no calls	# incorrect calls	Correct call rate %
	2C9*2	2C9*3	VKORC1						
1	wt/wt	wt/wt	GG	15	15	15	0	0	100
2	wt/wt	wt/*3	GG	15	14	14	1	0	93.3
3	wt/wt	wt/wt	GA	15	15	15	0	0	100
4	wt/wt	wt/wt	GG	15	15	15	0	0	100
5	wt/wt	wt/wt	GG	15	15	15	0	0	100
6	wt/wt	wt/wt	GA	15	14	14	1	0	93.3
7	wt/wt	wt/wt	GA	15	13	13	2	0	86.7
8	wt/*2	wt/wt	GA	15	14	14	1	0	93.3

9	wt/wt	wt/wt	AA	15	15	15	0	0	100
10	wt/wt	wt/wt	GA	15	15	15	0	0	100
11	wt/wt	wt/wt	GG	15	15	15	0	0	93.3
12	wt/wt	wt/wt	GA	15	15	15	0	0	100
13	wt/wt	wt/wt	AA	15	15	15	0	0	100
14	wt/wt	wt/wt	AA	15	15	15	0	0	100
15	wt/wt	wt/wt	AA	15	15	15	0	0	100
16	wt/wt	wt/wt	AA	15	15	15	0	0	100
17	wt/*2	wt/wt	AA	15	15	15	0	0	100
18	wt/*2	wt/wt	GA	15	15	15	0	0	100
19	wt/wt	wt/wt	AA	15	15	15	0	0	100
20	wt/wt	wt/wt	AA	15	15	15	0	0	100
21	wt/wt	wt/wt	AA	15	15	15	0	0	100
22	wt/*2	wt/wt	GA	15	14	14	1	0	93.3
cctl-1	wt/wt	wt/wt	GG	15	15	15	0	0	100
cctl-2	wt/wt	wt/wt	GA	15	15	15	0	0	100
cctl-3	wt/wt	wt/wt	AA	15	15	15	0	0	100
cctl-4	wt/wt	wt/wt	GG	15	15	15	0	0	100
cctl-5	wt/*2	wt/wt	GG	15	15	15	0	0	100
cctl-6	*2/*2	wt/wt	GA	15	15	15	0	0	100
cctl-7	wt/wt	wt/wt	GG	15	15	15	0	0	100
cctl-8	wt/wt	wt/*3	GG	15	15	15	0	0	100
cctl-9	wt/wt	*3/*3	GA	15	15	15	0	0	100
Total				465	459	459	6	0	98.7

Reproducibility per genotype:

Genotype	# samples tested [‡]	# correct calls	# incorrect calls	# no calls	Correct call rate %
2C9*2	465	461	0	4	99.1
2C9*3	465	464	0	1	99.7
VKORC1	465	464	0	1	99.7

[‡] Total tests for each genotype: 31 samples x 3 sites x 5 days repeats = 465 tests

Reproducibility per site per *Operator*:

Site/Operator	Locus	# Samples	# correct calls	# incorrect calls	# no calls	Correct call rate %
Site 1/Operator 1	2C9*2	93	91	0	2	97.85
	2C9*3	93	93	0	0	100
	VKORC1	93	93	0	0	100
Site 1/Operator 2	2C9*2	62	62	0	0	100
	2C9*3	62	62	0	0	100
	VKORC1	62	62	0	0	100
Site 2/Operator 3	2C9*2	93	93	0	0	100
	2C9*3	93	92	0	1	98.9
	VKORC1	93	92	0	1	98.9
Site 2/Operator 4	2C9*2	62	60	0	2	96.8
	2C9*3	62	62	0	0	100
	VKORC1	62	62	0	0	100
Site 3/Operator 5	2C9*2	155	155	0	0	100
	2C9*3	155	155	0	0	100
	VKORC1	155	155	0	0	100

Lot to lot reproducibility results

A study was conducted at the sponsor's facility to determine the reproducibility between kit reagent lots of the eQ-PCR™ LC Warfarin Genotyping assay. A panel of 24 whole blood samples was tested using three reagent lots for five days. The results per genotype for each lot are illustrated in the table below.

Loci	Lot	# samples tested	# tests [§]	Genotyping calls made	# correct calls	# incorrect calls	# no calls	Correct call rate %
2C9 *2	Lot 1	24	120	120	120	0	0	100%
	Lot 2	24	120	120	120	0	0	100%
	Lot 3	24	120	120	120	0	0	100%
2C9 *3	Lot 1	24	120	120	120	0	0	100%
	Lot 2	24	120	120	120	0	0	100%
	Lot 3	24	120	120	120	0	0	100%
VKORC1	Lot 1	24	120	120	120	0	0	100%
	Lot 2	24	120	120	120	0	0	100%
	Lot 3	24	120	120	120	0	0	100%

[§]Total tests for each Lot: 24 samples x 5 repeats = 120 tests

b. Linearity/assay reportable range:

Not Applicable

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

The sponsor conducted stability studies for DNA stored at -20°C and 2-8°C using 9 commercially available genomic DNA samples with 3 lots of the reagent kit. Stability studies for DNA stored at -20°C were conducted up to one year at 3-month intervals and for DNA stored at 2-8°C for 30 days. Results demonstrated consistent accurate results. The sponsor is continuing the stability studies for DNA stored at -20°C studies for 2 yrs. However, the sponsor recommends using DNA from fresh whole blood samples for genotyping assay.

Reagents for the Warfarin Genotyping Kit are shipped to the user laboratory on dry ice. The sponsor recommends storing all reagents at -20°C in dark. The kit is stable at -20°C for one year. After first use, all reagents should be stored at 2-8°C in dark. The sponsor states that under these conditions, reagents are suitable for use up to three months.

Controls: The sponsor does not provide positive or negative controls to be used with this assay. However, the sponsor has recommended that a positive (heterozygous and/or homozygous for the three genotypes) sample for each variant, a negative control (a sample that does not contain the mutation of interest, i.e., a wild type sample); and a “Non-Template Control” be included with each test run. All quality control requirements and testing should be performed in conformance with local, state and/or federal regulations.

The sponsor conducted preliminary studies to establish the Melting Temperature (T_m) cut-off values to identify the presence of individual SNP positions by the eQ-PCR™ LC Warfarin Genotyping Kit. A total of 64 individual samples were tested to evaluate the melting temperature Cut-Off value for each SNP type. The average melting temperature for each SNP was calculated and the Cut-Off value was determined based on Average value ± 3 SD. The SNP result for each sample was confirmed by bi-directional sequencing. Results are listed in the table below.

T_m values obtained for each genotype:

	VKORC1		2C9*2		2C9*3	
	Tm Low (variant)	Tm High (Wt)	Tm Low (variant)	Tm High (Wt)	Tm Low (Wt)	Tm High (variant)
Min:	53.31	62.36	54.52	63.47	53.26	60.32
Max:	55.92	64.52	58.44	66.74	56.49	63.16
Avg.	54.94	63.57	57.34	65.58	55.23	61.89
SD	0.80	0.65	0.90	0.59	0.90	0.86

The established Cut-Off value for each SNP:

SNP ID	Melting temperature for Wild type	Melting temperature for variant
VKORC1	63.57 ± 1.95oC	54.94 ± 2.40oC
2C9*2	65.58 ± 1.77oC	57.34 ± 2.70oC
2C9*3	55.23 ± 2.70oC	61.89 ± 2.58oC

d. Limit of Detection:

A study was conducted at the sponsor's facility to validate the recommended amounts of DNA required for genotyping test using eQ-PCR™ LC Warfarin Genotyping test. Six human genomic DNA samples at 7 different concentrations to obtain 0.1, 1, 10, 50, 100, 200 and 600 ng in the final reaction tube were used and 26 repeats conducted for each test level (total of 3312 tests). These six genomic DNA samples covered all 9 possible 2C9*2, 2C9*3, and VKORC-1 genotypes. Two of the six human genomic DNA samples were tested in two additional concentrations (300 ng & 400 ng) and repeated three times. Based on the results, there were two “no calls” at 0.1 ng DNA. The results of DNA limitation study are given in the tables below. For eQ-PCR™ LC Warfarin Genotyping assay, the sponsor recommends using a total of 50 – 200 ng genomic DNA in the final reaction tube.

Call rate based on DNA amount:

Input DNA Amount (ng)	Genotype	# samples	# repeats	# tests	# genotype calls	# correct calls	# incorrect calls	# no calls	Correct call rate (%)
0.1	2C9*2	6	26	156	468	466	0	2	99.6
	2C9*3	6	26	156					
	VKORC1	6	26	156					
1	2C9*2	6	26	156	468	468	0	0	100
	2C9*3	6	26	156					
	VKORC1	6	26	156					
10	2C9*2	6	26	156	468	468	0	0	100
	2C9*3	6	26	156					
	VKORC1	6	26	156					
50 [£]	2C9*2	6	26	156	468	468	0	0	100
	2C9*3	6	26	156					
	VKORC1	6	26	156					
100 [£]	2C9*2	6	26	156	468	468	0	0	100
	2C9*3	6	26	156					
	VKORC1	6	26	156					
200 [£]	2C9*2	6	26	156	468	468	0	0	100
	2C9*3	6	26	156					
	VKORC1	6	26	156					
300	2C9*2	2	3	6	18	18	0	0	100
	2C9*3	2	3	6					
	VKORC1	2	3	6					

400	2C9*2	2	3	6	18	18	0	0	100
	2C9*3	2	3	6					
	VKORC1	2	3	6					
600	2C9*2	6	26	156	468	468	0	0	100
	2C9*3	6	26	156					
	VKORC1	6	26	156					

[‡] The sponsor recommended DNA amount for eQ-PCR Warfarin Genotyping assay

Call rate based on genotype:

Allele	Genotype	# samples	# tests	# genotype calls	# correct calls	# incorrect calls	# no calls	Correct call rate (%)
2C9*2	wt/wt	4	740	740	739	0	1	99.9
	wt/*2	1	182	182	182	0	0	100
	*2/*2	1	182	182	182	0	0	100
2C9*3	wt/wt	4	728	728	728	0	0	100
	wt/*3	1	188	188	188	0	0	100
	*3/*3	1	188	188	188	0	0	100
VKORC1	G/G	3	552	552	552	0	0	100
	G/A	2	370	370	370	0	0	100
	A/A	1	182	182	181	0	1	99.5
Total			3312	3312	3310	0	2	99.9

e. Analytical specificity:

The sponsor conducted studies to demonstrate the interference of albumin (5 mg/ml), bilirubin (10 uM), triglycerides (4 mg/ml), heparin (14 U/ml), hemoglobin (80 mg/ml), magnesium chloride (6 mM), and lithium chloride (9 mM) and excess of extraction buffers on eQ-PCR assay. Nine blood samples covering all of the genotypes under evaluation were included. An aliquot of each sample was spiked with the above interfering substances and the results compared against the sample aliquot with no substance. The tests were repeated 6 times for each sample. The test results indicated no interference from above substances at the concentration levels tested and call rates remained unaffected.

f. Assay cut-off:

Not Applicable.

2. Comparison studies:

a. Method comparison with predicate device:

The sponsor conducted initial method comparison studies at 3 sites to demonstrate the accuracy of the device in comparison to bi-directional sequencing. This study included whole blood samples collected in EDTA from 150 donors. DNA was extracted at the sponsor's facility and aliquots of 50 unique DNA samples distributed to 3 sites (two external and the sponsor's facility) that participated in the study. The studies on the 150

samples resulted in 3 no-calls and 2 incorrect calls. The sponsor's study using 150 samples did not cover *3/*3 variant genotypes tested by the device. Therefore, the sponsor conducted a second study that included 2 additional *3/*3 samples. Genotyping results confirmed *3/*3 variant genotype for one of the samples, while the other new sample gave an incorrect genotype call of wt/*2,*3/*3 for the 2C9 gene. The sponsor performed a third study that included 7 new samples (4 *3/*3 and 3 *2/*2 samples). There were no incorrect calls or no-calls reported in this third study.

Genotype call-rates for three studies combined are given in table below.

Genotype	# of samples	# of correct calls	# of incorrect calls	# of no calls	Correct call rate
VKORC1	159	156	0	3 [‡]	98.11%
2C9*2	159	155	1 [‡]	3 [‡]	97.48%
2C9*3	159	155	1 [‡]	3 [‡]	97.48%

[‡] The incorrect genotype call reported a genotype of wt/*2, *3/*3.

[†] No-calls resulted from a system failure.

The genotype call-rate for each specific genotype based on all three studies combined is given in the table below.

Allele	Genotype	# samples	# genotype calls	# correct calls	# incorrect calls	# no calls	Correct call rate (%)	95%One-sided Confidence Lower Limit*
2C9*2	wt/wt	126	126	123	0	3	97.6	94.40%
	wt /*2	28	28	27	1 [‡]	0	96.4	87.68%
	*2/*2	5	5	5	0	0	100	47.98%
	Sub-total	159	159	155	0	3	97.5	94.59%
2C9*3	wt/wt	140	140	138	0	2	98.6	96.09%
	wt /*3	13	13	12	0	1	92.3	75.32%
	*3/*3	6	6	5	1 [‡]	0	83.3	54.28%
	Sub-total	159	159	155	0	3	97.5	96.40%
VKORC1	G/G	79	79	77	0	2	97.5	93.16%
	G/A	63	63	62	0	1	98.4	94.32%
	A/A	17	17	17	0	0	100	80.52%
	Sub-total	159	159	156	0	3	98.1	95.54%
Total for Assay			477	466	2 [‡]	9 [†]	97.7	96.18%

[‡] The incorrect genotype call reported a genotype of wt/*2, *3/*3.

[†] No-calls resulted from a system failure.

b. Matrix comparison:

The sponsor recommends only fresh whole blood samples only for use with this device. Therefore, Matrix comparison studies were not applicable.

3. Clinical studies:

a. Clinical Sensitivity:

Not Applicable.

b. Clinical specificity:

Not Applicable.

c. Other clinical supportive data (when a. and b. are not applicable):

4. Clinical cut-off:

Not Applicable.

5. Expected values/Reference range:

The following table provides the allele frequency across ethnic groups for the CYP2C9 and VKORC1 variants.

Expected Genotypic Frequency of Different Ethnic Groups:

	Caucasian ^{†1}	African ^{†2,3}	Asian ^{†1}
2C9*2			
wt/wt	78.7%	87.0%	100%
wt/*2	20.4%	8.7%	0%
*2/*2	0.9%	0%	0%
2C9*3			
wt/wt	88.0%	N/A	92.7%
wt/*3	11.6%	4.3%	7.3%
*3/*3	0.4%	0%	0%
VKORC1			
GG	39.1%	79.2%	1.9%
GA	46.7%	20.8%	18.3%
AA	14.2%	0%	79.8%

[†] The data presented are based on the genotypic frequencies in the following 3 publications:

1. Yuan et al., A novel functional VKORC1 promoter polymorphism is associated with inter-individual and inter-ethnic differences in warfarin sensitivity. Human Molecular Genetics 14(13): 1745-1751, 2005.
2. Kirchheiner and Brockmoller. Clinical Consequences of Cytochrome P450 2C9 Polymorphisms. Clinical Pharmacology & Therapeutics 77(1):1-16, 2005.
3. Schelleman et al., Warfarin response and Vitamin K Epoxide Reductase Complex in African Americans and Caucasians. Nature 415(6879): 428-432, 2007

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.